

# 07\_02 Cyclic Synthesis

MITTWOCH, 13.10.2021


## Goal-Setting

- Immobilized TdT-reaction with biotin-labelled primer and streptavidin-tagged magnet beads to create transitions

## Terms / abbreviations

- dNTP = Deoxynucleoside triphosphate

## Risk areas



## Required materials and / or information

- Chemicals:
  - 1x Buffer BXT Strep-Tactin XT elution buffer, IBA Life Science
  - 20 nM biotin-labelled primer, Ella Biotech
  - Endotoxin free water, Invitrogen
  - dNTPs, ThermoFisher
- Material:
  - 2 mL tubes (autoclaved), Eppendorf
  - 1000  $\mu$ L pipette tip
  - Lighter
  - Neodym Stabmagnet 4x60 mm, Maqna
  - PureCube HiCap StrepTactin MagBeads, Cube Biotech
  - Scissors

## Templates, devices, software

- None

## Preliminary work

- [07\\_01 Biotin Streptavidin Immobilization](#)
- Create a mastermix according to [03\\_01 ThermoFisher Protocol for TdT Tailing Reaction](#)

## Operation

- The following workflow is an example of a cyclic synthesis, it is not yet optimized
1. Add 100  $\mu$ L of the biotin-labelled primer and add 10  $\mu$ L of the magnet beads. Incubate for 10 min
  2. Add the magnet with the envelope and incubate for 5 min

3. Put the magnet with the envelope in a tube with endotoxin free water as a washing step for 5 min
4. Put the magnet with the envelope in a tube, which containing a mastermix for TdT reaction for 15 min
5. Washing step for 5 min (like in step 3)
6. Repeat steps 4 and 5 until the desired transitions are created
7. Dissolve the biotin-streptavidin bond by carefully pipetting 30  $\mu$ L of elution buffer around the envelope, stir, and incubate 10 min

## Troubleshooting

- Try to perform every experiment isolated from other experiments involving magnets, as magnet sticks can move and contaminate other experiments

## Disposal

- Autoclave trash bags, discard in S1 waste

## Follow-up work

- [📄 03\\_04 PCR for ssDNA Samples](#)