

# 00\_02\_10x-TBE-buffer-recipe

---

FREITAG, 14.5.2021

## Goal-Setting

- Manufacture of 10x Tris-borate-EDTA-buffer

## Terms / abbreviations

- EDTA = Ethylenediaminetetraacetic acid
- HCl = Hydrochloric acid
- NaOH = Sodium hydroxide
- PAA = Polyacrylamide
- TBE = Tris-borate-EDTA
- TRIS = Tris(hydroxymethyl)aminomethane

## Risk areas

 Hazard symbols



## Required materials and / or information

- Chemicals:
  - Boric acid, Merck
  - EDTA Lot, Sigma-Aldrich
  - MilliQ water, Sartorius arium pro VF
  - Tris Base, AppliChem
- Materials:
  - Beaker/Schottflask
  - Stirring fish
- All chemicals are solid substances
- Use TBE buffer only two times for gel electrophoresis

## Templates, devices, software

- Analysis balance, Kern ABJ 220-4NM
- pH meter, Knick Digital pH-Meter 646

## Preliminary work

- None

## Operation

- pH meter (ready to use in the chemical room)
  - If it is used for the first time, ask an experienced person around how to use it
  - Regularly calibrate the pH meter according to manufacturer
  - User guide:

- pH meter should measure about 6.5 when the solution is stored (small Falcon Tube)
- Take it out and unlock it (small switch at the top)
- Wash it carefully with some MilliQ water
- Hold into the solution and measure
- Wash again after use, turn off the switch and put it back into the small Falcon Tube

**The following recipe is standardized to 1 L**

1. Weigh the following components into a beaker or schottflask --> Tip: use at least a 1500 mL beaker (easier usage/measurement)
  - Tris Base: 121,14 g/L
  - Boric acid: 52,54 g/L
  - EDTA: 2,92 g/L
2. Fill up to 1 L with MilliQ water
3. Mix by stirring/shaking/inverting
4. Optimal pH = 8.3
  - a. If pH is too low, adjust with NaOH
  - b. If the pH is too high, do not adjust it with HCl, as this would cause errors with the PAA gel, because the Cl<sup>-</sup> ions will migrate first

## Disposal

- Observe all federal, state and local environmental regulations
- Here: Can be discarded in sink with a lot of water

## Troubleshooting

- None

## Follow-up work

- If a chemical runs out or is empty, please fill out the order list on the fridge next to the window
- [01\\_01\\_Agarose-gel-preparation](#)
- [01\\_03\\_Performing-agarose-gel-electrophoresis](#)
- [02\\_01\\_PAA-gel-preparation \(native\)](#)
- [02\\_01\\_PAA-gel-preparation \(denaturing\)](#)