

# 05\_05\_Colony-PCR

---

MITTWOCH, 2.6.2021

## Goal setting

- Extract and amplify the fragment from the picked colonies as preparation for Sanger sequencing

## Terms / abbreviations

- cPCR = Colony polymerase chain reaction
- PCR = Polymerase chain reaction

## Risk areas

- None

## Required materials/ information

- Chemicals:
  - 2x Taq PCR Mastermix, Qiagen
  - Forward Primer, Ella Biotech (M13 fwd 10  $\mu$ M)
  - Nuclease free water, ThermoFisher
  - Reverse Primer, Ella Biotech (M13 rev 10  $\mu$ M)
- Materials:
  - Ice
  - PCR tubes, Sarstedt
  - Pipettes, Eppendorf

## Templates, devices, software

- Nanodrop spectrophotometer, ThermoScientific Nanodrop 2000
- Thermocycler, Eppendorf

## Preliminary work

- [05\\_04\\_pGEM-T-Easy-Vector-System](#)

## Operation

### Detecting the right template DNA amount:

- General guidelines for low complexity DNA (e.g. plasmid, lambda or BAC DNA) are: 1 pg–10 ng per 50  $\mu$ L reaction volume
- If cDNA synthesis reaction mixture is used as a source of template, the volume of the template should not exceed 10% of the PCR reaction volume.
- To detect the desired amount, measure concentration in ng/ $\mu$ L according to [04\\_01\\_Spectralphotometer-Nanodrop](#)
- Then calculate how much is needed

### Prepare a Mastermix for the cPCR for needed amount of reactions:

- Always prepare at least for 1-5 reactions more than needed!

#### Composition of Mastermix



	A	B	C
1	Component	Volume for x reactions [μL]	x
2	2x Taq PCR Mastermix	500	50
3	M13 fwd 10 μM	100	
4	M13 rev 10 μM	100	
5	Nuclease free water	300	
6	Total	1000	

- Fill 20 μL of the Mastermix into the PCR Tubes (on ice!)
- Pick the colonies from the plates, streak them on a numbered master plate and put the pipetting tips into the PCR Tubes (on ice!)
- Shake the PCR Tubes carefully and remove the pipette tips

#### Run the following cPCR program:

#### cPCR Program



	A	B	C
1	°C	Time	Cycles
2	95	2'	
3	95	30"	30x
4	55	30"	
5	72	45"	
6	72	5'	

## Troubleshooting

- None

## Follow-up work

- Purify fragments with a cleanup-kit