

# 00\_01\_10%-APS-recipe

---

FREITAG, 14.5.2021

## Goal-Setting

- Recipe for 10% ammonium persulfate used for polyacrylamide gels

## Terms / abbreviations

- APS = Ammonium persulfate

## Risk areas



## Required materials and / or information

- Chemicals:
  - APS (solid), Roth
  - Autoclaved MilliQ water, Sartorius arium pro VF
- Materials:
  - 1.5 mL tubes, Eppendorf
  - Falcon Tube 50 mL, Sarstedt

## Templates, devices, software

- Analysis balance, Kern ABJ 220-4NM
- Vortexer, Scientific Industries (Vortex Genie 2)

## Preliminary work

- Search for gloves and safety goggles

## Operation

1. Weigh 10% (w/v) of APS into e.g. a Falcon Tube (e.g. 1 g APS for 10 mL H<sub>2</sub>O)
2. Add some autoclaved MilliQ water and fully dissolve the APS by vortexing
3. Fill with remaining volume of MilliQ water to reach the final desired volume
4. Mix by stirring/shaking/inverting
5. Aliquot 1-1.5 mL into Eppendorf Tubes and store at -20 °C in the iGEM freezer

## Disposal

- None

## Troubleshooting

- APS should always be stored in fridge
- For long-term storage, keep it in -20 °C freezer

## Follow-up work

- [02\\_01\\_PAA-gel-preparation \(native\)](#)
- [02\\_01\\_PAA-gel-preparation \(denaturing\)](#)