

# 03\_04\_Phusion-PCR-reaction

---

MITTWOCH, 2.6.2021

## Goal setting

- Performing a PCR reaction using the Phusion™ High-Fidelity DNA Polymerase

## Terms / abbreviations

- DMSO = Dimethyl sulfoxide
- PCR = Polymerase chain reaction

## Risk areas

- If spilled, always wipe surface with alcohol

### Hazard symbols



## Required materials/ information

- Chemicals:
  - 5x Phusion HF buffer, ThermoFisher
  - 5x Phusion GC buffer, ThermoFisher
  - 10 mM dNTP, ThermoFisher
  - 10 µM Forward primer, Ella Biotech
  - 10 µM Reverse primer, Ella Biotech
  - 50 mM MgCl<sub>2</sub> solution, ThermoFisher
  - Autoclaved MilliQ water, Sartorius arium pro VF
  - DMSO, ThermoFisher
  - Phusion DNA Polymerase (2 U/µL), ThermoFisher
  - Template DNA
- Material:
  - Ice or cooling rack
  - PCR tubes, Sarstedt
  - Trash bags, Th. Geyer GmbH & Co. KG

## Templates, devices, software

- Nanodrop spectrophotometer, ThermoFisher
- Thermocycler, Eppendorf
- Pipettes, Eppendorf

## Preliminary work

- Any reaction creating DNA that should be amplified
- [04\\_01\\_Spectralphotometer-Nanodrop](#) to detected desired amount of template DNA

## Operation

- Put the PCR tube on ice or the cooling rack
- Add the components in the order listed in the following table

Composition of Master Mix							
	A	B	C	D	E	F	G
1	Component	20 [μL]	50 [μL]	Final concentration	Unit	Linear PCR	Exponential PCR
2	MilliQ water	add to 20 μL	add to 50 μL			yes	yes
3	5X Phusion HF Buffer	4	10	1	x	yes	yes
4	10 mM dNTPs	0.4	1	200	μM	yes	yes
5	Forward primer 10 μM	1	2.5	0.5	μM	yes	yes
6	Reverse primer 10 μM	1	2.5	0.5	μM	no	yes
7	Template DNA *	x	x	0.2		yes	yes
8	Phusion High-Fidelity DNA Polymerase	0.2	0.5	0.02	U/μl	yes	yes

\* about 1 pg - 10 ng per 50 μl

- Run the PCR program in the thermocycler according to the following table

PCR Program for Thermocycler				
	A	B	C	D
1	Temperature [°C]	Time	Unit	Repeats
2	98	30	s	1x
3	98	5	s	30x
4	63	10	s	
5	72	30	s	
6	72	5	min	1x

If the PCR worked: increase 30 s to 1 min in step 5

## Disposal

- Autoclave trash bags, discard in S1 waste

## Troubleshooting

- Wear gloves to reduce the risk of DNase and RNase contamination
- The Phusion DNA Polymerase should be pipetted carefully and gently as the high glycerol content (50%) in the storage buffer may otherwise lead to pipetting errors

## Sources

 MAN0012393\_Phusion\_HighFidelity\_DNAPolymerase\_UG.pdf

## Follow-up work

- Purify fragments with a cleanup-kit
- If materials are empty, care about new order